
SO-APP+

Latest Version: 001.00.00
Date of publication: 2026/10/06

DA	Brugervejledning	5
DE	Gebrauchsanweisung	23
EN	Instructions For Use	41
FI	Käyttöohjeet	59
FR	Notice d'utilisation	77
IT	Manuale di istruzioni	96
NL	Gebruiksaanwijzing	114
NO	Bruksanvisning	132
SV	Bruksanvisning	150

DA	Dansk	DE	Deutsch	EN	English	FI	Suomi
FR	Français	IT	Italiano	NL	Nederland	NO	Norsk
SV	Svenska						

DA	Symboler	DE	Symbole	EN	Symbols	FI	Symboli
FR	Symboles	IT	Simboli	NL	Symbolen	NO	Symboler
SV	Symboler						

	DA	Medicinsk udstyr	DE	Medizinisches Gerät
	EN	Medical device	FI	Lääkinnälliset laitteet
	FR	Dispositif médical	IT	Dispositivo medico
	NL	Medisch hulpmiddel	NO	Medisinsk utstyr
	SV	Medicinteknisk utrustning		

	DA	Producent	DE	Hersteller	EN	Manufacturer
	FI	Valmistaja	FR	Fabricant	IT	Fabbricante
	NL	Fabrikant	NO	Produsent	SV	Tillverkare

	DA	Fremstillingsdato	DE	Herstellungsdatum
	EN	Reference number	FI	Valmistuspäivä
	FR	Date de fabrication	IT	Data di fabbricazione
	NL	Fabricagedatum	NO	Produksjonsdato
	SV	Tillverkningsdatum		

	DA	Referencenummer	DE	Bestellnummer
	EN	Catalog number	FI	Luettelonumero
	FR	Référence catalogue	IT	Numero di catalogo
	NL	Catalogusnummer	NO	Katalognummer
	SV	Referensenummer		

	DA	Serienummer	DE	Seriennummer
	EN	Serial number	FI	Sarjanumero
	FR	Numéro de série	IT	Numero di serie
	NL	Serienummer	NO	Serienummer
	SV	Serienummer		

DA	Se brugsanvisningen	DE	Lesen Sie Gebrauchsanweisung
EN	Read the instructions for use	FI	Lue käyttöohje
FR	Consulter le manuel d'utilisation	IT	Leggi il manuale utente
NL	Lees de gebruiksaanwijzing	NO	Les bruksanvisningen
SV	Läs användarhandboken		

DA	CE-mærke	DE	CE-Kennzeichnung	EN	CE marking
FI	CE-merkki	FR	Marquage CE	IT	Marchio CE
NL	CE-markering	NO	CE-merking	SV	CE-märkning

(01)

DA	UDI-DI	DE	UDI-DI	EN	UDI-DI
FI	UDI-DI	FR	UDI-DI	IT	UDI-DI
NL	UDI-DI	NO	UDI-DI	SV	UDI-DI

Anvendelsesområde

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The screenshot shows a mobile application form for patient feedback. At the top, there is a section titled "Følelse" (Feeling) with three circular icons representing different emotions: a sad face (red), a neutral face (yellow), and a happy face (green). Below this, there are two input fields: "Titel" (Title) and "Beskrivelse" (Description). The form is dark-themed with green and yellow accents.

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Scope of application

This User Guide contains all the information required to use the SO-APP+ mobile app safely and effectively. You must read this User Guide carefully and in full before using this device.

You are advised to keep this User Guide for future reference. If any of the information in this User Guide is unclear or if you have any questions, please contact France Développement Électronique (contact details provided at the end of this User Guide).

Description of the app

Intended Use

The SO-APP+ medical device is a mobile application that can be connected to an infusion pump in the SO-CONNECT+ range (V1 and 100) to extract and display data about infusions. SO-APP+ is designed for use by healthcare professionals to monitor patient treatment adherence, and also offers patients a different perspective on their prescribed treatment.

Performance

N/A

SO-APP+ distribution

The SO-APP+ app is distributed via the Apple and Google app stores:

For iOS devices	 App Store
For Android devices	 Google Play

Training

Healthcare professionals, certified homecare providers and clinicians can train individuals (patients and caregivers) who want to use SO-APP+ at home.

Intended users

The intended users are:

- Healthcare professionals trained in the use of the SO-APP+ medical device.
- Patients who have been taught to use SO-APP+ by a healthcare professional.
- For patients who cannot manage their treatment independently, caregivers who have been trained in the use of SO-APP+ by a healthcare professional.

	Healthcare professionals	Patients	Caregivers
Gender	All genders		
Age	No restrictions	18 years old or older	18 years old or older
Occupation	Licensed healthcare professional	No restrictions	No restrictions
Device users	Yes	Yes, provided the healthcare professional considers that the patient has the necessary mental and physical capacity to do so. Otherwise, no (in which case a caregiver is required).	Yes
Knowledge of the condition	Yes	Yes	Yes

Indications

The SO-APP+ mobile app is designed to help healthcare professionals monitor treatment adherence among patients using infusion pumps in the SO-CONNECT+ range.

The app retrieves and displays data from these infusion pumps, generating a quick overview of immunoglobulin administration over a period of up to 3 months. The healthcare professional can then see whether the patient is adhering to their treatment.

Active ingredient doses are expressed in ml and mg to facilitate reading and ensure alignment with prescriptions which are typically expressed in mg.

SO-APP+ is an additional tool that supports healthcare professionals in their day-to-day monitoring of patients.

Contraindications

The SO-APP+ app should not be used by patients with:

- Severe cognitive impairment
- Psychotic symptoms
- Or by patients lacking the necessary maturity (children).

Where patients cannot use the devices independently (e.g. children), they must be assisted by a caregiver trained in the use of the devices.

Patients with serious vision impairments should not use SO-APP+ without the assistance of a trained person.

Adverse effects

To date, no adverse effects associated with the use of the device have been identified.

Comments

If you need help setting up or using SO-APP+, please contact France Développement Electronique (contact details provided at the end of this User Guide).

You can also contact France Développement Electronique to report any unexpected issues or events.

Any serious incidents should be reported to the manufacturer and the competent authority of the Member State in which the user or patient resides.

For any questions concerning product quality, please contact the manufacturer at qualite@fde.fr

If you have any questions and/or comments about how the application uses your data, please contact the manufacturer at rgpd@fde.fr.

The contact details for the relevant authorities are available on the European Commission's website: https://health.ec.europa.eu/medical-devices-sector/new-regulations/contacts_en#nca

Warning

N/A

SO-APP+ service life

SO-APP+ has a service life of two years from the date it is available in the stores.

SO-APP+ characteristics

Classification	SO-APP+ is designated as a Class I medical device in compliance with Regulation (EU) 2017/745, Annex VIII, Chapter III, Rule 13.
Locking of restricted functions	Password-protected
Compatibility	Android and iOS
Service life	2 years

General overview of SO-APP+

The SO-APP+ device is designed for use by healthcare workers to check patient treatment adherence.
For patients, SO-APP+ offers an alternative way to track their treatment and gives them a fresh perspective of their SO-CONNECT+ pump.

Data protection

SO-APP+ only collects data provided by an SO-CONNECT+ pump or entered manually by the user. This data is not transmitted to a service provider or to the app developer.
Data is stored securely and exclusively on the smartphone. At no point is it transmitted to a data server or hosting provider.
When the app is uninstalled, all data is permanently erased.
Users can delete their data at any time via the application settings accessible through the host operating system.

SO-APP+ icons

Icon	Description
	Navigates back to the app's home page.
	Displays the app's main menu, providing access to SO-APP+ features.
	Shares pump information with a designated contact person.
	Indicates that the selected item can be edited. Tap to edit the item.
	Navigates to the left.
	Navigates to the right.
	Navigates to the first infusion available in the calendar.
	Navigates to the last infusion available in the calendar.

	Displays the Patient Information page.
	Displays the page containing infusion data.
	Shows that an infusion has been completed and gives access to the details. Green or red depending on the status of the infusion.
	An alert was triggered during the infusion.
	Feature for monitoring patient adherence.
	Information is available about the page you are currently viewing.
	A history log is available for this feature.

The Patient and Infusions icons are highlighted in green if the page is currently displayed on the screen.

Installing SO-APP+


SO-APP+ is available in both the Google and Apple app stores. Simply install the app as you would any other app.

There are two ways to find the app:

- Either by scanning the QR code found on the website and in this User Guide,
- or by searching for it in the store using the following keywords: So-App+.

Using SO-APP+ for the first time

When the app is launched for the first time, you must grant it access to the device's Bluetooth system. If permission is refused, you will not be able to use SO-APP+.

Bluetooth permission request



Once this step has been completed, you must select a user profile:

- Healthcare worker: this profile requires a password. It grants access to all SO-APP+ features (data retrieval, data display, data sharing, drug management, adherence, pump management).
- Patient: this profile grants access to basic features only (data retrieval, data display, data sharing, drug and profile management).

For each profile, the user's first and last name are requested for the sole purpose of customising the app and its features.


SO-APP+ home screen and main menu

Once the profile has been selected, SO-APP+ opens the application's home screen. This allows you to search for a pump or access other features via the menu in the top left.


Bluetooth enabled / Bluetooth disabled

Tapping the Menu icon brings up the following list:

Healthcare worker menu

Patient menu

Main menu [Patient, Healthcare worker]

From the main menu, you can:


- Scan a pump [Patient, Healthcare worker]
- Import a pump [Healthcare worker]
- View monitored patients [Healthcare worker]
- Change the app language [Patient, Healthcare worker]
- Edit your profile [Patient, Healthcare worker]
- View regulatory information about SO-APP+ [Patient, Healthcare worker]

A patient can edit the following information in their profile; however, only the dose and interval fields will be transmitted to the healthcare worker when the data is exported.

Patient profile


Pump connection [Patient, Healthcare worker]

SO-APP+ can detect SO-CONNECT+ pumps in Bluetooth-pairing mode within a range of under 5 metres.


Activating the export function on SO-CONNECT+

Tapping the SO-APP+ search button launches a 10-second pump-detection cycle.

Searching for pumps

If the app only detects one pump, data download begins automatically.
If multiple SO-CONNECT+ pumps are detected, you must select which pump to analyse from the list provided.


2 pumps found

Data retrieval [Patient, Healthcare worker]

When SO-APP+ starts extracting data from the pump, a dedicated download screen is displayed for the duration of the data transfer.

Once the download is complete, SO-APP+ will display the Patient Information screen.


All the data has been downloaded.

Data display [Patient, Healthcare worker]

The data downloaded by SO-APP+ is then organised into two categories, each with its own dedicated space within the app.


The Patient Information page is displayed immediately after the data retrieval phase and can then be opened by selecting the Patient icon.

On this page, users can modify treatment details, such as frequency and drug, by selecting the pencil-shaped Edit icons.

Patient page


The Infusion Information page can be accessed via the Infusions icon. This page includes a calendar listing all the infusions administered, as well as a button for adding new events to the calendar. It also features an information panel displaying the total volumes injected in ml over the last 30 and 90 days, together with the average doses in mg for the same periods.

Infusion page

Tapping the Syringe icon displays details of the completed infusion. These details include information about the volume, flow rate and all events recorded during the infusion.

By tapping the More icon, users can add an event which they think is important for tracking their treatment. To add the event, the user needs to enter its date, title and relevant details. Once confirmed, the event will be displayed in the calendar.

Infusion details

Data sharing [Patient, Healthcare worker]

After the data has been retrieved from the SO-CONNECT+ pump, it can be shared with others via the smartphone's existing messaging applications.


The app will prompt you to enter the patient-reported experience, the name of the file and any comments.

Feeling
 Title
 Description

Overall impression of the patient

The app supports two different formats:

- PDF: SO-APP+ generates a PDF listing the infusions with their details.
- SO-APP+: SO-APP+ generates a file in a proprietary format, enabling another healthcare worker to read the pump data as though it had just been collected.

File format


Each file can also be saved on the phone instead of being transferred to a third party. By default, the name of the generated file has the following format:

SerialNumber- YearMonthDay- HoursMinutesSeconds

Example: 000436-20250919-145223

Drug management [Patient, Healthcare worker]


SO-APP+ now lets you select the drug administered during the infusions from a list generated by the app.

Drug selection


This selection is made either after the data has been downloaded or from the Patient Information screen.

SO-APP+ also allows you to track the history if there have been any changes in the drug administered during treatment. To apply the change, you need to enter the start date of the new treatment.

Medication start date

SO-APP+ displays the active ingredient history when you select the corresponding icon on the Edit page.

Drug history


If you realise that you have made a mistake, SO-APP+ allows you to insert a new active ingredient into the history.

For a new pump which already has an existing history, the drug entered after downloading is applied to all infusions. The user must enter the treatment history. This history is automatically included when generating the treatment tracking report.

Monitoring adherence [Healthcare worker]


SO-APP+ allows healthcare workers to analyse patient treatment adherence more quickly through trend analyses based on multiple criteria. Whenever data is downloaded, the adherence feature provides the healthcare worker with a trend analysis using visual indicators (colour-coded smileys) generated from the list of events reported by the pump.

The healthcare worker can edit the adherence result based on their own criteria and the events displayed on the Adherence page. The app will prompt them to select the new status and to leave a comment if needed.

Different pages of the adherence feature

The adherence result is included on the first page of the PDF report.
The adherence criteria used by SO-APP+ can be viewed by selecting the Information icon at the top right of the Events page.


Event	Medium occurrence	High occurrence
Oedema	≥ 1	> 1
Redness	≥ 4	
Other effect	≥ 1	> 1
Flow change	≥ 1	
No treatment	≥ 7	≥ 14
Volume difference	< 20	≥ 20
Occlusion	≥ 2	≥ 4
System default		≥ 1

Adherence criteria

The values associated with the “No treatment” criterion are based on the treatment duration entered on the Patient Information page. A high occurrence corresponds to twice the average occurrence. Alerts are triggered with a ± 1 -day tolerance in order to avoid being overly restrictive.

Pump management [Healthcare worker]

SO-APP+ improves pump management compared with SO-APP, and now offers integrated multi-pump/multi-patient capabilities. You can now monitor all your patients' pumps in the app, even if you are not connected to the pump.

List of pumps already downloaded in SO-APP+

Every pump for which data is imported via the BLE connection is automatically saved in the memory. This means that you can review all the data from one pump at a later date.
By default, each added pump is named using the "Patient SerialNumber" pair, but you can change the name of a pump by selecting the Edit icon.

Likewise, if you export the data in SA+ format, your colleagues can view the data from this pump and add the pump to their list if they choose.
A patient can also send you their data in this format. Once it is imported, you will have access to the adherence feature described earlier in the document.

Analysis report

As we saw with the sharing function, SO-APP+ can generate a report in PDF format on request. This report includes the following information:

- Infusion history:
 - Patient’s assessment of their treatment
 - Type of event with severity levels plus date
 - Summary table of infusions
 - Comments from the nurse and doctor
 - Name given to the report
 - Description added to the report
 - Summary table of infusions
- The result of adherence (only in the Healthcare worker profile)

Analyze

Name : Patient 000508

Serial No. : SO-CONNECT+ V2 N° 000508

Period : from 2026/05/27 to 2026/08/27

General analyze opinion :

High priority events	Date	Medium priority events	Date
No drip during 14 days	2026/08/27 00:00		

Infusion History

Infusion interval (days) : 7

General treatment opinion : 

Configuration (Last):

Volume : 0 ml
Flow rate : 300 ml/h
Duration : 01:26:33
Infusion sites : 0

Patient's comments	Nurse's comments	Doctor's comments


Synthesis of infusions

Date	Volume infused - programmed	Flow max - min	Nb sites	Drug name	Alarm activation	Alarm missing settings	Events
Saturday 27/11/1999	0 ml - 300 ml	2.72 ml/h - 2.72 ml/h	0	Fluonate 200 mg/ml	0	0	/
Saturday 17/06/2000	200 ml / 0 ml 40.000 g / 0.000 g	300 ml/h - 2.72 ml/h	0	Fluonate 200 mg/ml	0	1	/
Monday 19/06/2000	31.24 ml / 0 ml 6.252 g / 0.000 g	40 ml/h - 40 ml/h	0	Fluonate 200 mg/ml	1	0	/
Sunday 23/07/2000	300 ml / 0 ml 60.000 g / 0.000 g	300 ml/h - 40 ml/h	0	Fluonate 200 mg/ml	0	0	/

SO-APP+ error codes

If a problem arises, SO-APP+ will display an on-screen error page with a message describing what has gone wrong and a button to go back to the main page. No error codes are reported by the application.

Should you experience any issues, please contact your distributor or France Développement Electronique (see contact details at the end of this User Guide) for guidance on the steps to take.

Error page

SO-APP+ unavailable

SO-APP+ can be uninstalled like any other app. All data stored in the app is deleted during the uninstallation process.

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